

FotoSan

Light activated disinfection

*Antibacterial treatment
without medication*

FotoSan indications

- Endodontics
- Gingivitis
- Periodontitis
- Periimplantitis
- Pericoronitis
- Caries


Safe and
without
side effects

FotoSan kills bacteria within seconds!

Improve your success
rate within root canal
fillings:



Fill the root canal with FotoSanAgent.




Activate FotoSan 30 seconds.



Rinse and dry canal.


Reduce bleeding,
swelling and pocket
depth:



Fill the pocket with FotoSanAgent.



Activate FotoSan 10 seconds in a deep pocket...



and/or 10 seconds on the outside of the pocket.

Designed, developed and manufactured in Denmark by: